



Certificate of Conformance

Certificate Number: EXMJ2250X27

Dated Issued: 03/04/2019

The blades were produced in accordance with Exsurco current design specification using only approved manufacturing processes and conforms to all established quality inspection and test criteria.¹

Product Item Number:	X100593
Item Description:	Pack, Amalgatome SK head, 5 pc
Item Lot Number:	X-0000000000027
Expiration Date:	01/2021

The blades were manufactured using a validated sterilization process:

STERIS Gamma Process Run ID:	102851A		
Processing Run Date:	27-JAN-2019		
Minimum Specified Dose (kGy):	25.0	Minimum Delivered Dose (kGy):	29.2
Maximum Specified Dose (kGy):	45.0	Maximum Delivered Dose (kGy):	36.0

Exsurco Medical is an ISO 13485:2016 certified and an FDA registered company. Exsurco Medical is committed to a policy of **QUALITY WITHOUT COMPROMISE**.

Rushell Ruff
Quality Department

03-04-2019
Date

[Signature]
Quality Manager

03-04-2019
Date

¹ Manufacturing data filed within the Quality Management System and available onsite.